Oxidase Reagent Droppers

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Italiano: pagine 3
Español: páginas 4

INTENDED USE

Oxidase Reagent Droppers are used in the Kovacs oxidase test as a qualitative reaction in the identification of non-fermenters and miscellaneous gram-negative bacteria.^{1,2}

SUMMARY AND EXPLANATION

The oxidase test is based on the production of an enzyme called indophenol oxidase. This enzyme oxidizes a redox dye (present in the reagent) which results in a color change of yellow to dark purple.^{1,2}

PRINCIPLES OF THE PROCEDURE

Indophenol oxidase, in the presence of atmospheric oxygen, oxidizes the phenylenediamine oxidase reagent to form a dark-purple compound, indophenol.^{1,2}

REAGENTS

Oxidase Reagent Droppers contain 0.5 mL of a 1% aqueous solution of N,N,N',N'-tetramethyl-p-phenylenediamine dihydrochloride which has been formulated with agents to insure maximum stability.

Precautions: *in vitro* Diagnostic

Avoid contact with the skin. Rinse thoroughly with water if spilled.

Storage Instructions Store at controlled room temperature 15 - 30°C (59 - 86°F). Protect from light. The reagent droppers need no refrigeration.

Product Deterioration: This reagent is hermetically sealed in an ampule, which affords protection of the solution from chemical instability until expiration date. Do not use if reagent is other than light yellow.

PROCEDURE

Materials Provided: Oxidase Reagent Droppers.

Materials Required But Not Provided: Ancillary culture media, filter paper, quality control organisms and laboratory equipment as required for this procedure.

Directions for Use:

1. Hold reagent dropper upright and **POINT TIP AWAY FROM YOURSELF**. Grasp the middle with thumb and forefinger and squeeze gently to crush ampule inside the dropper. **Caution: Break ampule close to its center one time only. Do not manipulate dropper any further, as the plastic may puncture and injury may occur.**
2. Tap bottom on tabletop a few times. Then invert for convenient drop-by-drop dispensing of reagent.
3. Add a few drops of oxidase test reagent to a strip of filter paper (Whatman No. 1 or equivalent). Streak a loopful of bacteria on the reagent-saturated paper with a platinum loop or wooden applicator stick. Do not add excess reagent, as it may cause the reaction to fade on oxidase-positive organisms.
4. Positive reactions turn the bacteria violet to purple immediately or within 10 to 30 sec. Delayed reactions should be ignored.

LIMITATIONS OF THE PROCEDURE

One should only allow 10 - 30 sec for a positive reaction. Any delayed reaction should be considered negative. Wire loop containing iron may give a false-positive reaction.

Perform the oxidase test on gram-negative bacilli **only** on colonies from **non-selective** and/or **nondifferential** media to ensure valid results.

Colonies to be tested must be isolated from other colonies. The use of 24 h cultures is recommended. Do not use refrigerated cultures without allowing them to reach room temperature.

Viscid colonies may be negative due to poor penetration of the reagent.

Reactions from weak oxidase-positive organisms, e.g., *Pasteurella* species, may be inaccurate. Results inconsistent with other biochemical reactions or with the organisms should be repeated.

User Quality Control Guidelines:

It is recommended that positive (*Pseudomonas aeruginosa* ATCC® 27853) and negative (*Escherichia coli* ATCC 25922) controls be run simultaneously with the organism to be tested.

AVAILABILITY

Description

Oxidase Reagent Droppers, packaged 50 droppers/carton.

For specific catalog number information, visit our website <http://www.bd.com/microbiology>, or contact the nearest Becton Dickinson Microbiology Systems office.

REFERENCES

1. Kovacs, N. 1956. Identification of *Pseudomonas pyocyanea* by the oxidase reaction. *Nature* (London) 178: 703.
2. Shigerei, J. 1992. Test methods used in the identification of commonly isolated aerobic Gram-negative bacteria, p. 1.19.16. In H.D. Isenberg (ed.), *Clinical microbiology procedures handbook*, Vol. 1. American Society for Microbiology, Washington, D.C.

TECHNICAL INFORMATION: In the United States, telephone Technical Services, toll free (800) 638-8663.

©1999 Becton Dickinson and Company.
ATCC is a trademark of the American Type Culture Collection.

Becton Dickinson Microbiology Systems
Becton Dickinson and Company
7 Loveton Circle
Sparks, Maryland 21152 USA

Becton Dickinson France S.A.
11 rue Aristide Bergès
38800 Le Pont de Claix, France

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